



*Spiracle Technology, U.S.A*



## FERNO2/Spiracle CPR Mask Directions for Use

The FERNO2/Spiracle CPR Mask provides a physical barrier between the rescuer and victim, eliminating direct contact of the rescuer's lips with unknown subject. ;START mouth-to-mouth ventilation without delay.

The Foldable Mask promotes an airtight seal to the face allowing ventilation through both the mouth and nose simultaneously. Proper training in ventilation and CPR should be obtained before using this device. Practice using on manikins so that you are ready for any emergency that might occur.

**NOTE: In the event of an emergency, always contact 911 (USA) or the local emergency medical service as soon as possible.**

**NOTE: In the event that you do not have the CPR Mask, start mouth-to-mouth ventilation without delay.**

---

### USE ON ADULTS AND CHILDREN:

1. Remove the mask/valve assembly from case. Push out the soft dome with your fingers and ensure that the one way valve is in place.
2. Put on protective gloves (optional on some models).
3. Clear the airway of any debris.
4. Place the mask over the patient's mouth and nose. Use the rim of the mask to retract the lower lip so that the patient's mouth remains open under the mask. Use thumbs to hold the mask and your fingers to grasp the lower jaw in front of the earlobes. Pull upward to open the airway and seal the mask.
5. Blow slowly into the airway valve twice and watch as the patient's chest rises. If the patient's chest does not rise; reposition the head and chin until you successfully ventilate the patient.
6. Ventilate the patient every five seconds for an adult or every three seconds for a child. Continue rescue breathing until the patient responds or you are relieved by a higher medical authority.
7. For use on a child simply reverse mask.

**NOTE:** If the patient vomits; remove the Foldable Mask and clean both the patient's airway and the mask with your gloved fingers. Resume rescue breathing promptly.

### USE ON INFANTS:

The FERNO2/Spiracle CPR Mask may be reversed to better fit an infants face. Place the nose section of the mask under the infants chin. Seal the Foldable Mask as described under the Adult / Child section and ventilate at least every three seconds using the previously described technique.

### USE WITH OXYGEN(OPTIONAL):

The FERNO2/SPIRACLE CPR MASK -MASK may be used with supplemental oxygen. Attach the oxygen tubing to the optional oxygen valve and set the flow rate. A head strap is included on all models that have the optional oxygen valve. If the patient is spontaneously breathing, use the head strap to position the Foldable Mask on the patient's face. If the patient is NOT breathing, begin mouth to mask ventilation immediately. The use of supplemental oxygen to the FERNO2/SPIRACLE CPR MASK -MASK will increase the delivered oxygen to the patient. The rescuer will not be affected by the use of supplemental oxygen.

**USE WITH OTHER VENTILATION EQUIPMENT:**

The FERNO2/SPIRACLE CPR MASK -MASK may be used with the Bag-Valve Mask or oxygen powered ventilation devices. Simply detach the valve/filter assembly and attach the Foldable Mask.



**DO NOT** pasteurize, boil, steam or gas autoclave this device. The "mask" component can be cleaned with cold liquid disinfectant solution only for same patient! This is a "single patient", "single use" item. **DO NOT REUSE PRODUCT.** Reuse of the product can be hazardous to the user and/or patient. Hazards like, but not limited to, degradation of the product and/or cross-contamination.

**TECHNICAL SPECIFICATIONS:**

Temperature Limits:

- Operating: -1°F (-18°C) to 104°F (40°C)
- Storage: -40°F (-40°C) to 140°F (60°C)

Inspiratory Resistance:

2.5 cm/H<sub>2</sub>O at 50 lpm flowrate

Expiratory Resistance:

3.5 cm/H<sub>2</sub>O at 50 lpm flowrate

Mask in Carrying Case:

Dimensions: 5-1/4" X 4-1/4" X 1-3/4"  
Weight: Approximately 4.5 oz.

LATEX FREE AND DEHP FREE

Re-Ordering Information Refer to below list

FERNO2/SPIRACLE CPR MASK -MASK RE-USABLE CPR MASK (Non-Autoclavable - clean by liquid chemical sterilization)		
ITEM	CODE	DESCRIPTION
1	T50520	FERNO2/SPIRACLE CPR MASK: Mask with valve, O2 Port, Head strap in Red Case
2	50521	Replacement one-way valve (Single-Use only) Fits above item

**CAUTION:** Federal Law restricts this device to use by or on the order of a physician.



Manufactured For:

**Spiracle Technology**  
10601 Calle Lee, #190  
Los Alamitos, CA 90720 U.S.A  
Tel : 714 418-1091 Fax : 714 418-1095  
E-mail : [info@spiracle.com](mailto:info@spiracle.com)  
[www.spiracle.com](http://www.spiracle.com)



EUROPEAN CONTACT:

Ferno Washington Italia s.r.l.  
via Benedetto Zallone n.26  
40066, Pieve di Cento (Bo) - Italy  
P +39 051 6860028 / F +39 051 6861508  
[info@ferno.it](mailto:info@ferno.it)



FDA Establishment Registration No.: 2030745

**LOT** 112227

